## Experiment 3 Grow your own microbes

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- |                                                                            |                       |
|----------------------------------------------------------------------------|-----------------------|
| • to pour = verser                                                         | to watch = observer   |
| • to stir/to mix= mélanger                                                 | to note = noter       |
| • to clean the spoon... = nettoyer la cuillère                             | to check = vérifier   |
| • to cover= couvrir/recouvrir                                              | to wash = laver       |
| • to measure the temperature of the water= mesurer la température de l'eau | the bag = le sac      |
| • dirty/clean = sale/propre                                                |                       |
| • the warm dark area = un endroit noir et chaud                            |                       |
| • the slice of bread = la surface du pain                                  |                       |
| • to wipe= frotter                                                         | to examine = examiner |
| • Watch the colour. What colour is it ? It's .....                         |                       |
| • ...                                                                      |                       |

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

#### Materials

**I will need:**

- **Sliced bread** (the cheaper the better) 
- **Sterile water** (boil the water in advance to sterilise it then allow to cool) 
- **Transparent sandwich bags** 
- **Permanent marker pen** 
- **Agar plates** – use these instead of bread 
- **Cotton buds** – these are used to sample surfaces 
- **Sellotape** 


	.....
	.....
	.....
	.....
	.....
	.....
	.....

### Method 1 (basic)

1. Predict which areas will be the most dirty/clean. Write in your book...
2. Use water to wet the slice of bread.
3. Wipe bread across the surface. For hands, wipe hand across the bread.
4. Place bread in sandwich bag and label with pen.
5. As a “control”, wet some bread and place in a bag **without** wiping a surface.
6. Tie bags shut and place in warm dark area (e.g., near a radiator) for several days. Check bread each day to see if mould has grown. It should look like the picture below after a few days.
7. Examine the bread and record results in the table.
8. Discuss which areas were dirty and which were clean.

### (Method 2 (advanced))

1. Tip: write on the agar side of the plate when marking plates (not the lid).
2. Predict which areas will be the most dirty/clean.
3. Wet the cotton bud with the water and then wipe cotton bud across the surface being sampled. Twist whilst you wipe to make sure the whole cotton bud is covered with microbes.
4. Remove lid from agar plate. Wipe cotton bud across entire surface of the agar using a “zigzag” pattern.  

5. For fingers: use marker pen to divide plate into four sections. Place a finger on the surface of the agar and press down gently.  

- 6.
7. Seal plates with Sellotape and place in warm dark place for 24-48 hours. Observe plates and record the results.
8. DO NOT remove the Sellotape or open the plates. Bacteria/fungi are dangerous and can cause infections.)

### Results table

Area		Predicted number of microbes OR "rank" (1 = most microbes)	Actual number of microbes OR "rank"
Toilet			
Desk			
Book			
Floor			
Hands (before washing)			
Hands (after washing)			
....			